

BLOOD BANK REFRIGERATOR



REFRIGERATOR

Krew Instrument's Blood Bank Refrigerator is a reliable and efficient solution for storing blood and blood components. Designed to maintain a consistent temperature range of 2°C to 6°C, it ensures that blood products are kept at the optimal temperature for safe and effective storage. With a range of features, such as an alarm system and a digital temperature display, Krew Instrument's blood bank refrigerator is an ideal choice for any blood bank or medical facility in need of a high-quality storage solution.

KEY FEATURES

- Capacity of 60 bags to 875
- Microprocessor controlled temperature controller.
- Interior illumination for working
- 7 Days circular chart recorder
- CFC Free Compressor.
- Specially designed stainless steel
- Outer Glass
- Lockable door for added safety of blood bags /samples.
- 180sec. compressor· ON' delay timer to safeguard the compressor.
- Password protected key.
- Stainless Steel Scratch Resistance roll out type drawers.
- High and low Alarm in case the inside temperature overshoots the preset.

CONSTRUCTION FEATURES

- Inner stainless steel (304) 0.8mm thickness and outer mild steel, 1.0mm thickness powder coated.
- PUF Insulation of 50 to 60 mm thickness.

SAFETY STANDARDS

TEMPERATURE SPECIFICATION

- Temperature range - 2 °C to 6°C
- Temperature accuracy - $\pm 0.5^{\circ}\text{C}$
- Uniformity - $\pm 1.0^{\circ}\text{C}$

TEMPERATURE CONTROL

- Operating 2°C to 6°C.
- Motor and blower arrangement to have uniformity of temperature in loaded
- Inside See-through unbreakable acrylic door to prevent loss of temperature.

DOCUMENTS

- Calibration certificate with supporting document of NABL Certificate (Govt. of India).
- User manual
- Warranty card

AVAILABLE CAPACITY

OPTIONAL ACCESSORIES

- Stabilizer
- Data Logger
- Front Glass Door
- UPS / Battery backup
- Ink less Chart Recorder

Model No	Inner Size (H x W x D) in CMS	Volume (in litres)	Approx No of Bags	No. of Baskets
KBR - 10	60 X 42 X 40	100	50	2
KBR - 11	80 X 50 X 50	200	125	3
KBR - 12	100 X 58 X 50	300	200	4
KBR - 13	90 X 63 X 63	350	250	4
KBR - 14	95 X 70 X 70	450	325	5
KBR - 15	110 X 75 X 715	600	425	5



 **Krew Instruments Pvt Ltd.**

Registered Address: No. 4/1101, Niraj Park,Wayle Nagar, Gandhare Road, Kalyan(w) Mumbai - 421301, Maharashtra, India

Manufacturing Unit : Gala No. 004, B-Wing, Ground Floor Shree Ganesh Compound, Pawar Section Chikhloli MIDC Ambernath (E) - 421506, Maharashtra, India

 +91 9224433419 / +91 8591239083 |  krewinstruments@gmail.com



EN 12469

MSME
MICRO, SMALL & MEDIUM ENTERPRISES