

# DONOR COUCH



Krew Instrument's **Donor Couch** is specifically designed for blood collection in blood banks and donation camps, ensuring the comfort and safety of donors during the extraction process. Equipped with a powerful motor system, the couch allows for quick shifting to shock position and immediate emergency measures, if needed. Designed with simplicity in mind, the models are user-friendly for phlebotomists and blood bank technicians. Krew Instrument's Donor Couch is a reliable and essential tool for any blood bank or donation center.

## KEY FEATURES

- Sturdy frame made of mild steel powder coated and covered with ABS body.
- KDC-1 has 2 motors for isolated movements of the back rest and leg rest.
- The chairs are supplied with 3 layered stands for placing the BCM; Pressure monitoring apparatus and other accessories.
- The chair occupies 30% less space than a conventional hospital bed.
- The wide arm rests can be optimally set at the desired height and also can be moved to the side.
- This allows vascular access at a relaxed posture, fewer episodes of hypotension and muscle cramps in patients or donors.
- In case a vasovagal attack occur, the donor's head can be lowered immediately and legs can be
- Lifted above his heart level easily; allowing blood flow back to the brain and other vital organs.
- The motors are silent and controlled with remote for each movements.
- 4 no's Heavy duty castor wheels are provided for easy movement of the chair and are lockable.
- Power supply: 230VAC 50 Hz
- Power Consumption : 100W
- Lifting capacity 200 kgs
- Battery Backup (Only for KDC-2 Model)

## AVAILABLE CAPACITY

Model No.	KDC-1	KDC-2
No of Motor	2	3
Seat W x D (inches)	24.5" x 26"	24.5" x 26"
Leg Rest W x D (Inches)	24.5" x 22"	24.5" x 22"
Back Rest W x D (Inches)	24.5" x 30"	24.5" x 30"
Hand Rest W x D (Inches)	8" x 18"	8" x 18"
When Open	7 feet	7 feet
Height from the ground	25"	25"
Adjustable Height	NA	25"+ 6"
Battery Back up	NA	2-3 Hrs



Krew Instruments Pvt Ltd.

**Registered Address:** No. 4/1101, Niraj Park,Wayle Nagar, Gandhare Road, Kalyan(w) Mumbai - 421301, Maharashtra, India

**Manufacturing Unit :** Gala No. 004, B-Wing, Ground Floor Shree Ganesh Compound, Pawar Section Chikhholi MIDC  
Ambernath (E) - 421506, Maharashtra, India

+91 9224433419 / +91 8591239083 | krewinstruments@gmail.com

