

PLATELET INCUBATOR (22°C)



Krew Instruments is proud to offer a cutting-edge Platelet Incubator (22°C) that provides precise temperature control for the safe storage of platelets. Our Platelet Incubator is designed to meet the unique needs of blood banks and transfusion centers, ensuring that platelets are stored in optimal conditions for maximum effectiveness. With its advanced features, including digital temperature control and alarms for temperature deviations, our Platelet Incubator provides peace of mind to laboratory professionals.

KEY FEATURES

- Capacity of 125 Liters. and 250 Liters.
- Microprocessor controlled Temperature Controller.
- Interior illumination for working area.
- CFC free compressor.
- Toughened Glass Door to monitor Agitation of platelet bags.
- 7 Days Circular Chart Recorder.

TEMPERATURE CONTROL

- Operating temperature of 22°C ± 1°C
- Inside acrylic door to prevent loss of temperature.
- Also internal fan switches off when the door is opened to avoid temperature.

CONSTRUCTION FEATURES

- Inner stainless steel (304) 0.8mm thickness and outer mild steel, 1.0 mm thickness powder coated.

SAFETY STANDARDS

- Lockable door for added safety of samples & materials.
- 180 seconds compressor· ON' delay timer to safeguard the compressor.
- Adjustable stainless steel Trays.
- High and low Alarm in case the inside temperature overshoots the preset value.

DOCUMENTS

- Calibration certificate with supporting document of NABL Certificate (Govt. of India).
- User manual
- Warranty card

OPTIONAL ACCESSORIES

- Stabilizer
- Data Logger
- UPS / Battery backup
- Ink less Chart Recorder

AVAILABLE CAPACITY

Model No	Size (H x W x D) in CMS	Capacity (in litres)
KPI - 125	46 x 59 x 46	125
KPI - 250	92 x 59 x 46	250



📍 Krew Instruments Pvt Ltd.

Registered Address: No. 4/1101, Niraj Park,Wayle Nagar, Gandhare Road, Kalyan(w) Mumbai - 421301, Maharashtra, India

Manufacturing Unit : Gala No. 004, B-Wing, Ground Floor Shree Ganesh Compound, Pawar Section Chikhholi MIDC Ambernath (E) - 421506, Maharashtra, India

☎ +91 9224433419 / +91 8591239083 | ✉ krewinstruments@gmail.com



EN 12469

MSME
MICRO, SMALL & MEDIUM ENTERPRISES