

WEIGHING SCALE



Krew Instrument's Micro Controller based Blood Weighing Scale accurately measures weight and volume of blood and its components, with the added convenience of easy conversion between the two. Its compare mode helps balance bags for refrigerated centrifugation, and it includes a zero set provision to indicate bag weight. Krew Instrument's Blood Weighing Scale provide accurate measurement of blood bags in grams and milliliters. The Weighing Scale comes with Abs moulded body with Tare button.

KEY FEATURES

- An LCD display of weight and volume with the accuracy of 1 gm./ml
- Can helps with a better balancing of the refrigerated centrifuge
- Easy prompt conversion of weight to volume
- Tare provision to all accounting for the weight of a blood bag
- An ABS moulded body
- Displays the volume and weight of all blood components
- Fine Auto Calibration
- Over load indication Capacity

TECHNICAL SPECIFICATION

- Input Voltage : 230 VAC
- Power Supply : SMPS
- Transducer : Load cell
- Weighing Range : 1 kg
- Accuracy : 1gm. / 1ml
- Calibration : Auto
- Display : LCD
- Weight : 1.5 kg

BLOOD BAG WEIGHING SCALE MODELS

Model No	SINGLE PAN WEIGH SCALE KWS – E	DUAL PAN WEIGH SCALE KWS – 2E
Description	<ul style="list-style-type: none"> • Accurate measurement of blood bags in Grams and Millilitre • Abs moulded body. • Tare Button. • Multikeys also available 	<ul style="list-style-type: none"> • Dual Pans for simultaneous measurement of 2 bags. • Accurate measurement of blood bags in Grams and Milliliters • Abs molded body. • Tare Button. • Components measurement on a single button press • Audio visual Alarm when the bags are balanced on both pans



 Krew Instruments Pvt Ltd.

Registered Address: No. 4/1101, Niraj Park,Wayle Nagar, Gandhare Road, Kalyan(w) Mumbai - 421301, Maharashtra, India

Manufacturing Unit : Gala No. 004, B-Wing, Ground Floor Shree Ganesh Compound, Pawar Section Chikhlohi MIDC
Ambernath (E) - 421506, Maharashtra, India

 +91 9224433419 / +91 8591239083 |  krewinstruments@gmail.com



EN 12469

संकेत : MSME
MICRO, SMALL & MEDIUM ENTERPRISES