

PHOTO STABILITY CHAMBER



Krew Instrument is a leading name in the pharmaceutical industry, manufacturing photo Stability Chambers as per ICH, WHO, and I.S standards. These chambers offer exceptional accuracy and reliability in controlling temperature and humidity, making them suitable for research applications that require a high level of control in a small footprint. Krew's photo Stability Chambers are specialized environmental test chambers that can achieve extreme conditions of temperature and humidity, adhering to strict guidelines and ensuring unparalleled precision in experimental conditions.

TECHNICAL SPECIFICATIONS

- Stand by cooling system: 1 unit
- Stand by Boiler system: 1 unit
- Control system: PLC-based system with touch screen HMI with PT 100 sensor: 1 unit
- 21 CFR Part 11 Compliant Desktop-Based Software: 1 unit
- Temperature + Humidity Scanner 4+4 channels with PLC/HMI Controller: 1 unit
- SMS Alert: 1 unit

MATERIAL OF CONSTRUCTION

- MS (mild steel)



Registered address: Unit No. 004, B-Wing, Ground Floor, Shree Ganesh Compound, Pawar Section, Chikhholi MIDC, Ambarnath (East), Maharashtra 421506

Call: +91 922 443 3419 / +91 859 123 9083 | **website:** [www,krewinstruments.com](http://www.krewinstruments.com)